wallace®

Sure-Pro® & Sure-Pro Ultra® Embryo Replacement Catheters

en

INSTRUCTIONS FOR USE

fr

MODE D'EMPLOI

de

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it

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es

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pt

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nl

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sv

BRUKSANVISNING

da

BRUGSANVISNING

no

BRUKSANVISNING

fi

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el

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Phone: +45 46 79 02 00

Fax: +45 46 79 03 02

CooperSurgical, Inc.
95 Corporate Drive
Trumbull, CT 06611 USA

CE2797

EC REP

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Sure-Pro® & Sure-Pro Ultra® Wallace® Embryo Replacement Catheters

INSTRUCTIONS FOR USE

These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient and/or clinician.

NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.

Wallace® embryo replacement catheters are single-use devices provided sterile for the introduction of embryos into the uterine cavity following *in vitro* fertilization.

WARNING:

- **READ THE COMPLETE DIRECTIONS FOR USE BEFORE USE.** Failure to properly follow the instructions, warnings and cautions may lead to serious surgical consequences or injury to the patient.
- **These procedures should only be performed by persons having adequate training and familiarity with these techniques.** Consult medical literature regarding techniques, complications and hazards prior to performance of these procedures.
- **TO BE USED BY, OR UNDER THE DIRECTION OF, QUALIFIED PERSONS IN LINE WITH LOCAL GUIDELINES GOVERNING *IN VITRO* FERTILIZATION, IF APPLICABLE.**
- **Discard if product or packaging is damaged.**

1 Description

These instructions are intended for use with sterile Wallace® embryo replacement catheters, product codes:

Code:	Description:	Inner Catheter type:
PP623	<i>Sure-Pro® Wallace® Embryo Replacement Catheter</i>	Classic, Soft, 23cm
PPB623	<i>Sure-Pro® Wallace® Embryo Replacement Catheter with Obturator</i>	Classic, Soft, 23cm
PPB623X	<i>Sure-Pro® Wallace® Embryo Replacement Catheter with Obturator (Inner catheter and outer sheath packed in individual sterile pouches)</i>	Classic, Soft, 23cm
PPS623	<i>Sure-Pro® Wallace® Embryo Replacement Catheter with Stylet</i>	Classic, Soft, 23cm
PPS623X	<i>Sure-Pro® Wallace® Embryo Replacement Catheter with Stylet (Inner catheter and outer sheath packed in individual sterile pouches)</i>	Classic, Soft, 23cm
PE623	<i>Sure-Pro Ultra® Wallace® Embryo Replacement Catheter</i>	Sure View**, Soft, 23cm
PEB623	<i>Sure-Pro Ultra® Wallace® Embryo Replacement Catheter with Obturator</i>	Sure View**, Soft, 23cm
PEB623X	<i>Sure-Pro Ultra® Wallace® Embryo Replacement Catheter with Obturator (Inner catheter and outer sheath packed in individual sterile pouches)</i>	Sure View**, Soft, 23cm
PES623	<i>Sure-Pro Ultra® Wallace® Embryo Replacement Catheter with Stylet</i>	Sure View**, Soft, 23cm
PES623X	<i>Sure-Pro Ultra® Wallace® Embryo Replacement Catheter with Stylet (Inner catheter and outer sheath packed in individual sterile pouches)</i>	Sure View**, Soft, 23cm

*Sure-Pro Ultra® with "Sure View**" inner catheter – for enhanced visibility under ultrasound.

The devices consist of a partially supported, flexible, end opening inner catheter and a detachable outer sheath that is attached to the inner catheter by a Luer compatible adapter.

The inner catheter is 16G with a lumen that is uniform throughout its length. It has a stainless steel support with a series of 1cm graduations at the proximal end and the flexible inner catheter protrudes from this by 8.5cm.

The outer sheath has a preformed angle, a series of 1cm graduations at the distal end and an adjustable marker to indicate orientation and depth of insertion.

The inner catheter protrudes from the outer sheath by 5cm.

The material of the Sure View® inner catheter is modified to contain small bubbles that improve ultrasound image quality.

The obturator and stylet provide differing levels of rigidity, and when attached to the outer sheath, in place of the inner catheter, form a smooth radius at the tip.

2 Indications for use

Introduction of embryos into the uterine cavity.

3 Contraindications

The catheter should not be used:

- In the presence of chronic cervical infection
- In the presence of or after recent pelvic inflammatory disease
- For intra-fallopian procedures

4 Cautions

- Low resolution ultrasound equipment, or use of 'Classic' embryo catheters may compromise the visibility of the catheter.
- When performing transabdominal ultrasound, ensure that the patient has a full bladder.
- The patient should be fully prepared for embryo transfer before preparation and loading of the catheter commences.

5 Instructions for Use

- 5.1** Place patient in the lithotomy, dorsal, or left lateral position. The cervix should be exposed with a speculum and gently swabbed with cotton wool moistened with normal saline or medium.
- 5.2** When using transabdominal ultrasound follow normal procedures such that an optimal image of the uterine cavity is achieved.
- 5.3** The adjustable marker on the outer sheath can be set to indicate the desired cervical length/depth. Care should be taken when moving this marker to ensure that the orientation notch still corresponds with the direction of the pre-formed angle on the outer sheath.
- 5.4 For two stage transfer technique only –**
lock the hubs of the obturator/stylet and the outer sheath together and advance through the cervical canal to the inner os. If slight resistance is experienced, withdraw the catheter. Using thumb and forefinger curve the stylet/obturator to complement the patients anatomy and pass the catheter in this position to the internal os. Gently uncouple and remove the obturator/stylet leaving the outer sheath in place.

WARNING

The stylet/obturator and outer sheath hubs must be locked together. Removing the stylet/obturator from the outer sheath, prior to placement, may result in damage to the endometrial tissue and bleeding.

- 5.5** Under sterile conditions, fill a 1ml syringe with medium and eject any air bubbles. Remove the tip protector and firmly attach the syringe to the Luer compatible adapter of the inner catheter. Flush the catheter through with medium, remove syringe, and refill with medium. Reconnect the syringe and expel medium until the 0.1ml mark is reached, with the catheter tip pointing upwards to facilitate the removal of air bubbles.
- 5.6** Place the end of the catheter carefully into the medium droplet containing the embryos and inject a small amount of medium to break the boundary of surface tension that may appear.
- 5.7** Aspirate the embryos into the catheter and hand to person performing embryo transfer.
- 5.8a For single stage transfer technique only –**
with the inner and outer hubs locked in position, advance the catheter so that the inner catheter passes through the external and internal os, into the mid-uterine cavity.
- 5.8b For two stage transfer technique only –**
the inner catheter should be advanced through the outer sheath already positioned (**5.4 above**) so that the inner catheter passes through the external and internal os, into the mid-uterine cavity.
- 5.9** It may be necessary to twist the catheter as it negotiates the cervical canal. When using ultrasound manipulate the probe to obtain the best image of the inner catheter.

WARNING

The catheter should never be forced against digitally felt resistance while inside the uterine cavity, as forcing the catheter may result in damage to the endometrial tissue and bleeding.

The outer sheath, stylet or obturator should not be advanced further than the internal os, and should certainly never enter the uterine cavity, as this may result in damage to the endometrial tissue and bleeding.

The inner catheter should not be pulled back further than the last graduation. Removing it from the outer sheath completely may result in damage to the endometrial tissue and bleeding.

The supported inner catheter is only to be used with the Sure-Pro® outer sheath as provided. The supported inner catheter should not be used without the outer sheath as this may lead to injury to the patient.

- 5.10 For single stage transfer technique only –**
If slight resistance is experienced, withdraw the catheter and advance the outer sheath until only the tip of the inner catheter is visible. Aligning the most distal marking on the base of the inner catheter with the hub of the outer sheath ensures a smooth radius protrudes. Whilst held in this position, pass the sheath and catheter to the internal os and then advance the inner catheter through the internal os, into the mid-uterine cavity. In the rare event where the cervical canal cannot be negotiated by this method, stylets are available to assist insertion. Use code No. PPS623/PPS623X or PES623/PES623X and then follow the two stage transfer technique.
- CAUTION: The Sure-Pro® product range is not compatible with product code 1816NST (the outer sheath provided with product code 1816NST will not facilitate use of the supported inner catheter).**
- 5.11** Place tip of catheter approximately 1cm from the fundus. Twist the catheter to remove any mucus from the tip.
- 5.12** Expel the embryos.
- 5.13** Leave time for the uterine environment to equilibrate and slowly remove the catheter.
- 5.14** Hand catheter back to embryologist to check that no embryos remain. Remove the speculum and allow the patient to rest.
- 5.15** Dispose of the catheter in accordance with local medical hazardous waste practices.

The above are only suggested instructions for use. The final decision on the technique used is the responsibility of the clinician in charge.

Sure-Pro®, Sure-Pro Ultra®, Sure View® and Wallace® are registered trademarks of CooperSurgical, Inc. and its subsidiaries.
CooperSurgical is a registered trademark of CooperSurgical, Inc.
© 2019 CooperSurgical, Inc. All rights reserved.

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en

Single use.

Do Not Reuse: Medical devices require specific material characteristics to perform as intended. These characteristics have been verified for single use only. Any attempt to re-process the device for subsequent re-use may adversely affect the integrity of the device or lead to deterioration in performance.

Sterile unless unit container is opened or damaged.

Destroy after single use. Do not resterilize.

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