

# New IVF Workstations

For the Proper Environmental Maintenance of Living Cells.





## The New IVF Workstations

- Vertical Laminar Flow (available in 3, 4, 5, and 6-foot lengths)
- Integrated Warming Plate, with Customized Heating Patterns
- External Water Circulator & Heater
- Stereo Scope Installation
- Less than 55db for Decreased Noise Distraction
- VOC Air Filtration (optional)
- Electrical Receptacle Integrated to Back Plane
- Floor Stands (available in 30" & 36" heights) & Adjustable Stand
- Expert On-site Installation
- One-year Warranty (includes all parts and labor)

For ordering information Call your MidAtlantic Devices Representative at  
**Call 1-800-648-1151**

## Workstation Accessories

- Heating Patterns for Warming Plate Choose from the popular T-85L and T-45L warming plate heating patterns, (shown below) or from ten additional patterns or customize your IVF Workstation to your precise usage requirements.



- T-85L (Left) Shown  
T-85R (Right)
- T-45L (Left) Shown  
T-45R (Right)
- Custom Pattern

### • Video Systems, Displays, VCR's & Printers



- Flat Screen Mounting:  
Flush, Articulating Arm,  
Inside, Outside and  
Custom Positions



- Cameras
- Digital, CCD Systems
- Custom Installations



- VOC Filtration System
- Flexible
- Secure
- Ergonomic



- Water Heater/Circulator
- Quiet
- Reliable
- Visual Monitoring

### • CO2 Manifold/Humidifier Systems



- Single Flow Meter  
Gas flow meter 0-29  
liters/hour with  
regulator valve. Single  
gas manifold  
(not shown).



- Dual Flow Meter  
2 flow meters 0-29  
liters/hour. Dual  
regulator valves.  
Dual gas manifold  
(not shown).

- Triple Flow Meter  
Available



- Mini-Incubators  
Integrated to IVF  
Workstations  
or Free Standing.



- Accessories For optimizing temperature, ph and humidity environments. Available as glass hoods and aluminum blocks to fit any IVF dish or tube.