

Micropipet Type:

- ICSI
- Holding
- Assisted Hatching
- PZD
- Denuding
- Biopsy
- Stem Cell
- Piezo Drill
- Custom

Customization

We offer customization of our standard products so you may slightly modify any existing micropipet design. For example, you may want a longer spike on an ICSI micropipet, or a longer distance from the pipet tip to the placement of the bend. We will gladly make micropipets to your specifications. When placing your order, please let us know your preferences.

If you would like to completely design a new micropipet, you must complete this form and then fax it to us at 001 (434) 295-5912 or to your distributor. Custom micropipets must be used for research only and are not intended for the IVF market.

Specifications:

- _____ Length of overall pipet (up to 9cm)
- _____ Length of taper (cm)
- _____ Inner diameter of opening (µm)
- _____ Bevel (Yes or No)
- If yes:
 - _____ Spike or
 - _____ Polished or
 - _____ Neither
- If no:
 - _____ Flat and polished
 - _____ Flat and not polished
- _____ Requested Angle, if any (15°, 20°, 25°, 30°, 35°, 40°, 45°)
- _____ Angle placement from tip
- _____ Sterile (Yes or No)

Please provide a drawing of your custom(ized) micropipet:

Custom features:

Intended Use:

Contact Information:

Name

Organization

Address

City

State

Zip

Country

Telephone

Fax

Email

Research Tool

These tools should be used according to their intended use in officially approved research projects.