

NEW

Glass Denudation System

Get the handle that fits from the First Name in Micropipets



Priced
Competitively

Silicon bulbs
included with
pipets

ERGONOMICALLY DESIGNED HANDLE

- Lightweight
- Comfort grip
- Auto-clavable for maximum sterility

BOROSILICATE GLASS PIPETS

- Exact smooth edge for gentle denudation
- Precise inner diameter
- Full range of sizes (μm)

Denuding Micropipets

Denuding micropipets are used to remove the cumulus cells from the oocyte. These micropipets are made from 1.5 mm O.D. glass, have a blunt opening and are sold straight. They are available in a variety of sizes. All pipets are CE marked, FDA cleared, sterilized by gamma radiation and MEA and LAL tested.

Denuding Micropipets (10 per box)

Reference No.	Description	I.D.
MDP-120-0	Straight Denuding Micropipet, with blunt opening	120 μ m +/- 5 μ m
MDP-130-0	Straight Denuding Micropipet, with blunt opening	130 μ m +/- 5 μ m
MDP-135-0	Straight Denuding Micropipet, with blunt opening	135 μ m +/- 5 μ m
MDP-140-0	Straight Denuding Micropipet, with blunt opening	140 μ m +/- 5 μ m
MDP-150-0	Straight Denuding Micropipet, with blunt opening	150 μ m +/- 5 μ m
MDP-160-0	Straight Denuding Micropipet, with blunt opening	160 μ m +/- 5 μ m
MDP-170-0	Straight Denuding Micropipet, with blunt opening	170 μ m +/- 5 μ m
MDP-175-0	Straight Denuding Micropipet, with blunt opening	175 μ m +/- 5 μ m
MDP-190-0	Straight Denuding Micropipet, with blunt opening	190 μ m +/- 5 μ m
MDP-250-0	Straight Denuding Micropipet, with blunt opening	250 μ m +/- 5 μ m

Denuding Handle

Reference No.	Description
HANDLE-BL	Black Denuding Handle

Same space-saving packaging used for all micropipets

