


IVF Pasteur Pipets

- 5.75" (14.6cm) and 9" (22.9cm) standard sizes
- Sterilized by dry heat
- Available with or without cotton plug
- Mouse embryo and LAL tested
- Extensively cleaned and baked to remove endotoxins
- Compact packaging allowing removal of one pouch at a time
- Made from Borosilicate glass

Specifically
Prepared for
the IVF Lab



IVF Pasteur Pipets

IVF Pasteur Pipets, specifically prepared for the IVF laboratory, are available in standard sizes of 5.75" (14.6cm) and 9" (22.9cm). Pipets can be purchased with or without cotton plugs. They are sold in boxes of 90 (30 pouches of 3 each) or cases of 1,000 (20 trays of 50 each). We have designed convenient packaging for the boxes of 90 to include a perforated opening on each box which allows you to pull a pouch out easily while leaving the remainder of the pipets protected. This compact, but sturdy, packaging means more space on the lab shelf.

Pasteur Pipets

| Reference No. | Description |
|----------------|---|
| PP-5.75-90 | Box of 90, 30 pouches of 3 pipets each, 5.75" (14.6cm) long |
| PP-9-90 | Box of 90, 30 pouches of 3 pipets each, 9" (22.9cm) long |
| PP-5.75-1000 | Case of 1000, 20 trays of 50 each, 5.75" (14.6cm) long |
| PP-9-1000 | Case of 1000, 20 trays of 50 each, 9" (22.9cm) long |
| PP-5.75-90PL | Box of 90, cotton plugged, 30 pouches of 3 pipets each, 5.75" (14.6cm) long |
| PP-9-90PL | Box of 90, cotton plugged, 30 pouches of 3 pipets each, 9" (22.9cm) long |
| PP-5.75-1000PL | Case of 1000, cotton plugged, 20 trays of 50 pipets each, 5.75" (14.6cm) long |
| PP-9-1000PL | Case of 1000, cotton plugged, 20 trays of 50 pipets each, 9" (22.9cm) long |